



2 POSSIBLE ways to obtain your report :



1

Automatically :

At the end of the test, click **request a report** and in 3 minutes it will be sent to you automatically.

2

Manually :

Identify patient data file : Find the correspondence ID-Name of the patient by launching your MindPulse test and clicking on the button **Open folder**

Each file is identified by: Year-month-day-hour minutes_first part mail clinicien_patient ID .bin.

Open your folder containing patient data files by clicking **View patients list**

Drag your patient data file into your email or copy and paste to your computer desktop so you can easily recall it and send it to us (by the email address you gave us as your username). Object **«Request for report»**. NB: if you send an email using an email address other than your username, please remember your identifying email address.

Send us this email with your patient's data file attached to:

 **resultats@mindpulse.net**

NB : On line You must forward the email with the data file received from your patient to this same address : **resultats@mindpulse.net**

Data files are encrypted and anonymized by a unique ID.

In order not to damage them NEVER TRY TO OPEN THEM and NEVER CHANGE THE NAME OF THE FILE