



Instructions to install MindPulse test MAC VERSION

Before installation, please
ensure that you have :

- A **connected 64 bit computer**.
- A good quality **mouse with a cord** (no cordless mice)
- One of these **operating MacOS versions** : High Sierra, Mojave, Catalina, Big Sur.

NB: standarized data was performed with a mouse with a cord. Depending on the availability of a mouse with a cord and the patient's ease, one can use the touchpad.

Before installation, we
encourage you to :

- Close other running software.
- Mute the sound.
- Adjust the brightness (high) while ensuring the quality of the contrasts and comfort of the patient.
- Ensure an optimal distance from the screen (50 to 70 cm) as well as having good posture.

Contents

I | RAPID INSTALLATION p. 2

II | ASSISTANCE IN INSTALLATION p. 3-7



I. SIMPLE INSTALLATION

1 By clicking on the Apple logo, **download** the installation zip file

2 **Double-click** on the downloaded installation file.

3 Start the installation by clicking and **maintaining it**

4 Choose the save **destination** and click.

5 Choose the installation **type** and click **install**.

6 Successful installation, **Mindpulse is ready !**



II. ASSISTANCE IN INSTALLATION

Alert and safety messages


Warning messages by your antivirus are simply related to the originality of MindPulse. It is not yet recognized by all operating systems or antivirus programs. As well, your computer may trigger a security alert : **We guarantee that our download platform www.mindpulse.net/ is completely secure.**

1

By clicking on Apple logo, **download** the installation zip file

2

Double-click on the downloaded installation file.




Mac warns you that it does not know the application.
Click **OK**.




3

Click **apple menu** located at the top left.




4

Choose **System Preferences**.



5


Click **Security & Privacy**.






6

Click **General** tab and **Open anyway**.



7


Click **Open**






8

Start installation by pressing **Continue**.



9


Choose the save **destination** and click **Continue**.





10

Choose the installation **type** and click **install**.



11

Successful installation, **Mindpulse is ready !**

